



Grass[®] EC3[®] Conductive Adhesive Gel

Accuracy, ease of use & patient comfort

EC3 Conductive Adhesive Gel is a smooth, translucent gel for attaching electrodes and conducting signals during EEG and Sleep studies. EC3 is easy to apply, gentle on the skin and simple to remove.

Overview:

- Water-soluble: easy clean-up with warm water; reduced potential to stain clothing
- Smooth consistency for easy application
- Sets quickly to maintain position
- Hypoallergenic

For optimal results during long-term EEG studies, cover each electrode with a small piece of flexible adhesive wrap.

natus[®]

Natus Medical Incorporated
3150 Pleasant View Road
Middleton, WI 53562 USA
Tel: 1-800-356-0007
1-608-829-8500
Fax: 1-608-829-8709
natus.com

Ordering Information

Grass EC3 Conductive Adhesive Gel

| Part Number | Type | Size | Quantity |
|-------------|-------------------------|-------------------|----------|
| EC3 | Conductive Adhesive Gel | 1.76oz (50g) tube | 12/box |

Learn more about the full line of Natus Neuro products by contacting your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**

Convenient online portal (US Customers only)

Natus Neuro Store – www.natusneurostore.com